Abb. 1: Y-Zugang zum distalen Radius von radiopalmar  
Fig. 1: Y-shaped skin incision over the radiopalmar aspect of the wrist


Abb. 2: Spalten des 1. Strecksehnenfaches  
Fig. 2: Opening of the first extensor compartment


Abb. 3: Spalten des 3. Strecksehnenfaches und Herausheben der langen Daumenstrecksehne  
Fig. 3: Opening of the third compartment and retraction of the extensor pollicis longus

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Fig. 5: The palmar ridge of the first extensor compartment must be removed to assure optimum surface contact of the plate to the bone.

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Fig. 6: The distal screw holes are drilled by using the fixed angle centering sleeve. The holes for the fixed angle connections are tilted in an angle of 10° distally to the axis of the flap. This tilt should be considered when applying the centering sleeves as well as during insertion of the fixed angle screws.

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Fig. 7: Distal positioning of the correction plate

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
Fig. 8a

Correction in sagittal plane  
Present situation:  $-22^\circ$  Aim:  $+10^\circ$   
→ Correction angle:  $32^\circ$   
→ Bisector:  $16^\circ$


Fig. 8b

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Fig. 8 a/b: Firstly a drawing is made in correspondence with the current x-rays for the correction osteotomy in the sagittal plane. The dorsal tilt of the articular surface of the radius (f.e.  $22^\circ$ ) is determined through a lateral x-ray and the comparison with the x-ray of the non-injured side, which shows the desired position of the articular surface (usually  $10^\circ$  palmar tilt). The planned angle of correction (f.e.  $32^\circ$ ) is defined as the sum of the dorsal tilt and the desired palmar tilt. The angle of the osteotomy (f.e.  $16^\circ$ ) is the bisector of the correction angle.

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Fig. 9 a/b: A drawing is made to define the angle of correction and the angle of the osteotomy for the ulnar inclination. The x-ray in the frontal plane of the healthy forearm shows the desired ulnar inclination (f.e.  $25^\circ$ ). The correction angle (f.e.  $9^\circ$ ) results from the difference of the desired ulnar inclination to the present ulnar inclination of the injured side (f.e.  $16^\circ$ ). The osteotomy angle (f.e.  $4,5^\circ$ ) is the bisector of the correction angle.

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Fig. 8 c/d: The angle formed by the distally fixed plate with the radial shaft corresponds to the malalignment of the articular surface of the radius. After opening osteotomy the proximal part of the plate adapts to the radial shaft when the distal fragment assumes its anatomical position.

Correction of ulnar inclination  
Present situation:  $+16^\circ$  Aim:  $+25^\circ$   
→ Correction angle:  $9^\circ$   
→ Bisector:  $4,5^\circ$


Fig. 9c

Fig. 9d

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Fig. 9 c/d: The height ulnarely of the bone graft is set by the ulna variance.

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Fig. 10: After marking the designated site for the osteotomy the plate is removed. According to the drawing the osteotomy is made with an oscillating saw. The osteotomy is made proximal of the fixed angle screw holes in the bisector of the intended angle of correction in both planes.

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Fig. 11: After the osteotomy is completed the plate is distally fixed to the radius. The osteotomy is opened up and therefore the proximal part of the plate is reset onto the shaft of the radius. Temporary fixation with plate holding forceps. Adjusting the length of the radius and fixation of the plate under x-ray control.

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Fig. 13: Finally the bone graft is fixed with a lag screw

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Fig. 14 a: Pre-operative situation with shortening and dorsal tilt of the radius. Clearly visible arrosion of the lunatum because of the traumatic ulna impaction syndrome.

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Fig. 14 b: Early post-operative situation. The interpositioned bone block can be seen well.

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Fig. 14 c: Final result (1 year post-operative). The radius is restored in form and length. Reversion of the changes on the lunatum.

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Fig. 15 a: Intraoperative situation after the corrective osteotomy of the radius. On special request of the patient the interpositioning of the bone graft was avoided.

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Fig. 15 b: The osteotomy gap is "empty" early post-operatively. The angle stable fixation allows an early mobilisation.

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Fig. 15 c: Shortly after the explantation of the osteosynthesis material. The radius is restored in form and length. The gap of the osteotomy is filled with new bone even without bone grafting.