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Fig. 1: The drill guide is assembled (26-350-03 into 26-350-01) and is then carefully positioned on the bone to be fixed.

ATTENTION:

It is important to carefully remove any soft tissue overlying the bone, since the points on the tip of the guide must be engaged in bone, both to prevent 'wandering' and to ensure accurate length calibration. The position and alignment of the guide may be checked on an image intensifier, as these will determine the insertion point and direction of the guide wire.

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Fig. 2: Once the drill guide is in optimum position, the reduction sleeve (26-350-05) is introduced into the drill guide in order to control accurately the position of the guide wire (26-351-00 or 26-350-00). The guide wire is advanced slowly, under x-ray control, to ensure that it is optimally positioned perpendicular to the fracture/fusion, in the mid-axis of the bone. Its tip should enter, but not penetrate the opposite cortex of the bone. The optional tip-threaded guide wire (26-350-00) may be used to reduce the risk of the wire backing out during subsequent stages of the procedure.

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Fig. 3: A 1 mm fixation wire (26-351-00) is now passed through one of the six peripheral holes of the collar of the guiding handle (26-350-01); this pin may be advanced across the fracture/fusion site under x-ray control, in order to maintain a good reduction during insertion of the screw. When its position is judged satisfactory, the protruding end of the fixation wire should be bent or cut short, so that it does not obstruct access to the guide. At the end of the procedure it is either removed, or else cut flush to the bone surface and left in position to provide two-point fixation.

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Fig. 4: Once the fixation wire is in the correct position, the measuring sleeve (26-350-07) is passed over the protruding end of the guide wire to determine the required drilling depth.

ATTENTION:

To ensure accurate length measurement, the tip of the drill guide must remain firmly applied to the surface of the bone, and the measuring sleeve must be pushed over the guide wire as far as possible.

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Fig. 5: The reduction sleeve is now removed from the drill guide and the cannulated drill (26-350-09) is passed over the guide wire, the drill stop having been set at the depth indicated by the measuring sleeve.

The drill is carefully advanced along the guide wire under x-ray control, until the drill stop reaches the end of the drill guide, indicating that the drill has reached its full depth.

ATTENTION:

It is recommended that the drill is withdrawn and cleaned one or more times during insertion, to prevent it from 'binding' to the guide wire; should this occur, the guide wire is likely to be withdrawn with the drill, in which case it should be repositioned (using the reduction sleeve if necessary) before moving to the next stage (Step 6).

Insertion of the guide wire and of the cannulated drill is best done using a suitable, slow-speed, small power instrument. If not available, then insertion of the cannulated drill may be done by hand, using the drill handle (26-350-11) contained in the set.

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Fig. 6: The screwdriver (26-350-17) is used to lift a screw (26-300-xx or 26-310-xx) of correct length from the rack in the tray.

The screw is passed over the guide wire and is inserted until increased resistance is felt as the trailing thread enters the bone

ATTENTION:

Should increased compression be required in order to produce optimum fixation (for example in the case of a nonunion) then the high compression screw (26-310-xx) may be selected in preference to the standard screw (26-300-xx).

Due to interfragmentary compression we recommend to select a screw length which is 2 mm shorter than the adjusted drill depth.

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Fig. 7: At this stage the guide wire should be removed, the screwdriver re-inserted and the screw advanced until the screwdriver abuts on the end of the drill guide.

ATTENTION:

The guide wire should always be removed before the screw is fully tightened; otherwise, its tip may become bent, making it impossible to remove.

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Fig. 8: The drill guide is now removed, and the screw is advanced a further one or two turns, in order to ensure that the trailing thread is well buried beneath the surface of the bone. If there is any difficulty in re-engaging the screwdriver, this may be facilitated by loosely repositioning the guide wire and using it to centre the end of the screwdriver in the socket of the screw.

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
Step 1: Drill guide

The drill guide is assembled (26-370-03 into 26-350-01) and is then carefully positioned on the bone to be fixed. The two notches on the drill guide enable either a closer or more distant positioning of the handle on the proximal fragment.

ATTENTION:

It is important to carefully remove any soft tissue overlying the bone, since the points on the tip of the guide should be engaged in bone to ensure accurate length calibration.

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
Step 2: Temporary fixation wire (optional)

A 1-mm fixation wire (26-351-00) may be passed through one of the six peripheral holes of the collar of the guiding handle (26-350-01); this pin is advanced across the fracture/fusion site under x-ray control, in order to maintain a good reduction during insertion of the screw. When its position is judged satisfactory, the protruding end of the fixation wire should be bent or cut short, so that it does not obstruct access to the guide. At the end of the procedure it is either removed, or else cut flush to the bone surface and left in position to provide 2-point fixation.

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

Step 3: Drilling

Once the fracture is accurately repositioned and the drill guide is in correct position, the drill (26-370-09) is inserted, either by hand or preferably, using a suitable power driver. Insertion of the drill under x-ray control (image intensifier) ensures precise positioning. Once the drill has reached the required depth, the calibration is read, in order to determine the correct length of the implant.

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Step 4: Insertion of the screw
The screwdriver (26-370-17) is used to lift a screw (26-320-xx) of correct length from the rack in the tray.
The screw is inserted until the screwdriver abuts on the end of the drill guide.

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Step 5: Optimum compression
The drill guide is now removed, and the screw is advanced a further one or two turns, in order to ensure that the trailing thread is well buried beneath the surface of the bone, whilst at the same time achieving optimum compression.