Operationstechnik – proximale Polfraktur des Skaphoids Surgical Technique – Proximal Pole Fracture of the Scaphoid

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Fig. 1: The drill guide is assembled (26-370-03 into 26-350-01) and is then carefully positioned on the bone to be fixed. The two notches on the drill guide enable either a closer or more distant positioning of the handle on the proximal fragment.

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ATTENTION:

It is important to carefully remove any soft tissue overlying the bone, since the points on the tip of the guide should be engaged in bone to ensure accurate length measurement.

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Fig. 2: Once the drill guide is in optimum position, the reduction sleeve (26-375-05) is introduced into the drill guide in order to control accurately the position of the guide wire (26-375-00). The guide wire (0.8 mm) is advanced slowly, under x-ray control, to ensure that it is optimally positioned perpendicular to the fracture/fusion, in the mid-axis of the bone.

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Fig. 3: Once the fixation wire is in the correct position, the measuring sleeve (26-350-07) is passed over the protruding end of the guide wire to determine the required drilling depth. The usual screw lengths are 18–22 mm.

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ATTENTION:

To ensure accurate length measurement, the tip of the drill guide must remain firmly applied to the surface of the bone, and the measuring sleeve must be pushed over the guide wire as far as possible.

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Fig. 4: The reduction sleeve is now removed from the drill guide and the cannulated drill (26-375-06) is passed over the guide wire, the drill stop having been set at the depth indicated by the measuring sleeve. The drill is carefully advanced along the guide wire under x-ray control, until the drill stop reaches the end of the drill guide.

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NOTE:

Insertion of the guide wire is best done using a suitable, slow-speed, mini power unit. During the drilling process, the drill should be moved back and forth several times, in order to avoid too much collection of saw dust inside the cannulation at the tip of the drill. The insertion of the cannulated drill may be done by hand, using the drill handle (26-350-11) contained in the set.

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Fig. 5: The screwdriver (26-375-17) is used to lift a screw (26-322-xx) of correct length from the rack in the tray. The screw is passed over the guide wire and is inserted.

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Fig. 6: At this stage the guide wire should be removed, the screwdriver re-inserted and the screw advanced until the screwdriver abuts on the end of the drill guide.

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ATTENTION:

The guide wire should always be removed before the screw is fully inserted; otherwise, its tip may become bent, making it impossible to remove.

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Fig. 7: The drill guide is now removed, and the screw is advanced a further one or two turns, in order to ensure that the trailing thread is well buried beneath the surface of the bone. If there is any difficulty in re-engaging the screwdriver, this may be facilitated by loosely repositioning the guide wire and using it to centre the end of the screwdriver in the socket of the screw.

Operationstechnik – Fingergelenksarthrodese Surgical Technique – DIP Arthrodesis

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Fig. 1: A central y-shaped incision is made over the distal end phalanx.

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Figs. 2 a/b: As a preparation for the following end-phalanx arthrodesis, a 0.8-mm guide wire (22-636-08) that has been pointed at both ends is first inserted under visual control from the distal interphalangeal joint into the distal phalanx. Following correct repositioning, the wire is then advanced into the middle phalanx from distal to proximal.

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Figs. 3 a/b: This procedure facilitates placing the wire in a central (mid-axis) position as required for correct screw insertion.

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Fig. 4: In order to determine the correct screw length the HBS Mini screw is placed intra operatively over the corresponding end phalanx. The leading thread should be positioned in the center of the middle phalanx. The usual screw lengths are about 20 mm.

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Fig. 5: The cannulated drill (26-375-06) is inserted into the distal cortex of the middle phalanx.

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Figs. 6 a/b: The cannulated HBS Mini screw (26-322-xx) is inserted from distal via the guide wire.

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Figs. 7 a/b: The small diameter of the screw guarantees stable fixation while preventing bone burst in the distal phalanx.

Ø 3,55 mm