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Fig. 1: Ulnoradial impingement syndrome with severe instability of the proximal ulnar stump following Sauve-Kapandji procedure after several previous surgical operations on the distal radioulnar joint.

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Fig. 2: Surgical exposure through the old scar, which is usually slightly extended in proximal direction, thereby identifying and preserving dorsal sensory branch of the ulnar nerve. This is followed by exposure of the fused ulnar head and the proximal stump between the extensor carpi ulnaris and extensor digiti quinti muscles. Longitudinal opening of a relatively dense scar between these two structures is required as well.

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Fig. 3: Spherical head reamers are now used to ream the proximal surface of the fused ulnar head.

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Fig. 4: The depth and shape of the reamed socket is checked with an appropriate spherical head gauge under fluoroscopic control.

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Fig. 5: The spherical socket should embrace two-thirds of the prosthetic head. However, the reamed area should be kept as shallow as possible to prevent arthrodesis fracture.

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Fig. 6: The medullary canal is opened with a pointed awl at the distal end of the ulnar stump, then a small rasp is hammered in as far as it goes. Depending on the planned size of the prosthesis, the appropriately sized rasp should be subsequently inserted.

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Fig. 7: Insertion of the trial prosthesis. Its conical end should be located slightly proximal to the fused ulnar head.

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Fig. 8: To prevent fractures in the thin bone edge of the ulnar head, we recommend indirect reduction of the new socket to the spherical head by performing an oblique, incomplete radial osteotomy.

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Fig. 9: The radius osteotomy is carried out proximal to the arthrodesis.

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Fig. 10 a/b: After widening the bone gap with a spreader clamp, the depth and shape of the socket are checked with an appropriate spherical head gauge with fluoroscopy. If the trial prosthesis has been inserted, a trial head can also be used for this purpose.

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Fig. 11: The definitive stem of the prosthesis is inserted with the conical drive-in tool.

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Fig. 12: To prevent fractures in the thin wall of the ulnar head socket, the prosthetic head is inserted, having the radial osteotomy distracted on the ulnar side in order to "distalize" the fused ulnar head.

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Fig. 13 a/b: As soon as the bone gap has been closed on the radius, the prosthetic head should perfectly articulate with the spherical joint socket (given correct preoperative length planning). A space of 1–2 mm will then form postoperatively between the prosthetic head and the joint socket.

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Fig. 14 a/b: The osteotomy is fixed with a dorsoulnar small-fragment plate. Local scar tissue can be used as a flap to cover the neck and head of the prosthesis. This flap can be sutured transosseously to the fused ulnar head. The tendon sheath of the extensor carpi ulnaris muscle is approximated to the sheath of the extensor digiti quinti muscle. This dorsalizes the extensor carpi ulnaris, thus providing additional prosthesis coverage.

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Follow-up Treatment

An ulnar Scotch Cast® U splint is worn until the wound has healed. Thereafter, a removable wrist splint is worn for another three weeks following suture removal.