Operationstechnik – Ulnaverkürzung Surgical Technique – Ulnar Shortening

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Fig. 1: Incision is made on the ulnopalmar side under tourniquet control with the forearm positioned dorsally. This ensures a good cover of soft tissue over the compression plate.

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Fig. 2: Access through the Septum intermusculare between the flexor and extensor muscles until the ulnar bone is exposed. Putting aside of the M. pronator quadratus in the distal part of the bone.

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
Fig. 3: Precise adaption of the compression plate to the palmar side of the ulnar bone is achieved by bending. This prevents tilting of the bone plate after the osteotomy.

Fixation on each side of the estimated position of the osteotomy, without tightening the screws completely. However it is important, that the screws are inserted at the proximal side of the compression hole, in order to allow for a maximum sliding distance and the rotational stability is secured.

In case of strong resistance of the screw during implantation, it can be necessary to tap even with self-tapping screws.

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Fig. 4: Loosening of the screws and lifting of the plate. Selection of the saw guide, which is depending on the width of the osteotomy (3 mm, 4 mm or 6 mm). The guide is then placed beneath the plate. One of the fixation notches is placed in the 3rd hole from distal and another in the distal part of the middle sliding hole.

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Fig. 5: The saw guide enables a parallel osteotomy directly underneath the 4th hole from distal.

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
Fig. 6: Fixation of the distal screws after removal of the bone disc.

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
Fig. 7: Closure of the osteotomy gap through traction of the
plate and possible Ulnar-duction. Tightening of the proximal
screws.

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
Fig. 8: The remaining proximal and distal holes are filled with
screws in an eccentric position for further compression. Finally
insertion of a lag screw in the middle hole.

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Fig. 9: Ulnar impaction with an ulnar + situation of 2 mm. Arthroscopical image of a broad tear of the TFCC.

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Fig. 10: Ulnar shortening of 4 mm with the special compression plate. X-ray control after 3 months which shows a united osteotomy. No clinical complaints.

In case of a **shortening of the radius** the plate is inserted in the same manner from the radiopalmar side. The osteotomy is carried out across the cancellous area. The insertion of a lag screw is not necessary. Therefore only a 6-hole plate is required.

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