
## Operationstechnik – Ulnaverkürzung Surgical Technique – Ulna-Shortening

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**Fig. 1:** Incision is made on the ulnopalmar side under tourniquet control with the forearm positioned dorsally. This ensures a good cover of soft tissue over the compression plate.

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
**Fig. 2:** Access through the Septum intermusculare between the flexor and extensor muscles until the ulnar bone is exposed. Putting aside of the M. pronator quadratus in the distal part of the bone.

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
**Fig. 3:** The shortening plate of the ulna (26-163-07) is placed palmarly at the inner edge of the ulna. From the distal end of the ulna the plate is positioned approx. 3 cm proximal. Should the plate not have optimal surface contact in this area, it is mandatory to adapt the plate precisely to the palmar surface of the ulna by bending. This prevents tilting of the bone after the osteotomy.

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**Fig. 4 a:** The connected drill guides allow for orthogonal bicortical drilling with a drill of a maximum diameter of 2.0 mm. They can be taken off after drilling with the tip of the screwdriver (26-950-17). The screws are inserted first of all **monocortical** in the sequence shown above. Please be aware, that the 3rd distal hole is also drilled but not filled with a screw at this stage. Alternatively either the non-locking 2.7 mm T-Drive screw (26-901-xx) or the locking 2.7 mm T-Drive screw (26-502-xx) can be used. Should you prefer the larger screws of 3.5 mm (26-903-xx or 26-504-xx) please note, that before drilling the drill guides need to be taken off. This is necessary, since a larger drill does not fit through the drill guides.

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**Fig. 4:** In order to determine the optimal screw length with the round depth gauge (25-219-04) please subtract 2 mm from the measured length. Should you use the flat depth gauge (26-945-00), which is included in the set, then please use the side of the depth gauge which shows the plate and add 2 mm.

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
**Fig. 5:** Loosening of the screws and lifting of the plate. Selection of the saw guide, which is depending on the width of the osteotomy (3 mm, 4 mm or 6 mm). The guide is then placed beneath the plate. One of the fixation notches is placed in the 3rd hole from distal and another in the distal part of the middle sliding hole.


By tightening of the screws the saw guide is fixed during the sawing.

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**Fig. 6:** The saw guide enables a parallel osteotomy directly underneath the 4th hole from distal. The saw blade should be not thicker than 0.5 mm, not wider than 15 mm and have a length of approx. 40 mm in order that the saw can swing freely in the slots. After sawing the screws can be loosened and the saw guide is taken out.

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**Fig. 7:** After removal of the parallel cortico-cancellous bone wedge the 3rd screw from distal is inserted and all distal screws are fixed **bicortical**.

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**Fig. 8:** Closure of the osteotomy gap through traction of the plate and possible Ulnaruktion.

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**Fig. 9:** The two proximal screws in the long holes are tightened bicortical, in order to fixate the closure of the osteotomy.

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**Fig. 10:** In the next step the middle proximal hole is drilled and filled. A lag screw is inserted in the remaining hole above the osteotomy. This leads to additional stability of the osteotomy.

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In case of a **shortening of the radius** the plate is inserted in the same manner from the radiopalmar side (alternat. a dorsal approach is possible). The osteotomy is carried out across the cancellous area. The insertion of a lag screw is not necessary. Therefore only a 6-hole plate (26-163-06) is required. The saw guide is not used with this plate type.